



Since 2021

SHIVANSH MEDICAL SERVICES

“YOUR HEALTH IS OUR WEALTH”

Regd. Address:- G 09 VARDHMAN MARKET SEC TOR-02 DWARKA NEW DELHI-110075

E-Mail : shivanshmedicalsolutions@gmail.com

Contact us : +91-8920692107 or +91-9350604227

DL No: 148031/148032

IE Lic.: BEAPS3883Q

GST No: 07BEAPS3883Q1ZO

ISO Certification: 9001:2015

MSME: UDYAM-DL-10-0039768

COMBAT KITS



IFAK KITS



ENHANCE MEDICAL KIT



PARAMEDIC KIT

Combat first aid kits are vital components of a soldier's gear, designed to provide critical medical support in the demanding and often unpredictable environments of combat and operational scenarios. Unlike conventional first aid kits, military versions are engineered with specific challenges of combat in mind. These kits are constructed to be lightweight, compact, and fit seamlessly in a soldier's loadout without adding unnecessary bulk. One of the key factors when it comes to military first aid kits is their adaptability. The contents address a variety of injuries and medical conditions: tourniquets, bandages, dressings, and personal protective equipment are some of the essential items included.

These kits can be customized based on the specific requirements of a mission, the environment in which soldiers will operate, and the level of medical training available.

Since its combat introduction by the U.S. Army in 2005, the Combat Application Tourniquet® - a true one-handed tourniquet - has decreased the mortality rate due to extreme exsanguination. The U.S. Army Institute of Research deemed it to be 100% effective.

Now, C·A·T Resources has developed the all-new GEN7. **CAT** Following extensive field reviews, the **CAT** GEN7 has been re-designed, making it the most effective Combat Application Tourniquet® to date. The new GEN7 is simpler to use thanks

to a single application protocol, single buckle routing system and further improvements to its components and construction.

GEN7



Identical models available in black [tactical], orange [First Responder], and blue [trainer].

WHEN EVERY SECOND COUNTS

Despite the proven effectiveness of the Combat Application Tourniquet®, the new **CAT** GEN7 goes even further toward making it easier and faster to use. All six components of the Combat Application Tourniquet® have been re-designed and improved on the **CAT** GEN7.

1. SINGLE ROUTING BUCKLE

Decreased blood loss, effective slack removal, fewer windlass turns, simplified training with single protocol application standards

2. WINDLASS ROD

Increased diameter, enhanced strength, aggressive ribbing improves grip

3. WINDLASS CLIP

Bilateral beveled entry, rapid windlass lock, bilateral buttress, added strength

4. WINDLASS STRAP

New reinforced strap, gray color for tactical considerations (Black only)

5. STABILIZATION PLATE

Reinforced, beveled contact bar, improved comfort

6. FREE-MOVING INTERNAL BAND

Patented band within band, truly even distribution of circumferential pressure



FAST. SAFE. EFFECTIVE. THE GEN7 ADVANTAGE

The Single Routing Buckle is a direct response to in-depth reviews of field and testing reports. Featuring a new single-routing buckle system, the **CAT** GEN7 allows for simpler and faster application with more efficient slack removal resulting in decreased blood loss with fewer windlass turns.



Manufactured by:
C·A·T Resources, LLC
483 Lakeshore Pkwy
Rock Hill, SC 29730
803.325.9300



HYFIN[®] VENT Chest Seal Twin Pack

The HyFin[®] Vent Chest Seal Twin Pack from North American Rescue sets the standard for the prevention, management and treatment of an open and/or tension pneumothorax potentially caused by a penetrating chest trauma. Two seals, each individually packaged, are designed to allow the user to apply a seal to both an entry or an exit wound while also giving the option to only apply one and store the other until needed as the situation requires. This innovative design provides 3-vented channels designed to prevent airflow into the chest cavity during inhalation while allowing air to escape through the vent channels during exhalation. The 3-vent channels are designed to allow blood to escape and also provide a backup fail-safe system, as even if two of the three channels become obstructed, the vent is designed to remain fully operational.

Advanced adhesive technology provides superior adhesion as the HyFin[®] Vent Chest Seal Twin Pack is specifically designed to overcome adverse conditions in austere environments where the casualty may be covered in blood, sweat, body hair, or other environmental contaminants when a reliable chest seal is mandatory. Packaged in a more compact, rugged, easy-to-open twin foil pouch, each chest seal also includes a gauze pad to wipe the wound surface prior to application. The HyFin[®] Vent Chest Seal Twin Pack has a large red pull tab for single-step peel-and-apply application, while the transparent backing facilitates easy placement over the wound area, conforming to the patient's chest.

NAR packages each chest seal using its signature Red Tip Technology[®] red tear notches in a rugged, easy-to-open, foil pouch to ensure high quality performance. The HyFin[®] Vent Chest Seal Twin Pack meets or exceeds the current EMS Standard of Care and TCCC & TECC Guidelines for treatment of penetrating injuries to the chest and is the superior prehospital chest seal.

Features:

- Designed for the prevention, management and treatment of an open and/or tension pneumothorax potentially caused by a penetrating chest trauma
- Patented, new design with 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of both entry/exit or multiple penetrating injuries
- Advanced adhesive technology designed to provide a superior seal in the most adverse conditions, including sweaty or hairy casualties
- Easy-to-grip, large red tab for single step, peel-and-apply application that allows for the burping of the wound if necessary

Specifications:

- Device: H 6 in. x W 6 in. (excluding tab)
- Packaged Folded: H 7.5 in. x W 4.5 in. x D 0.25 in.
- Packaged Unfolded: L 7.5 in. x W 9 in. x D 0.13 in.
- Weight: 2.5 oz.
- US Patents: 7,504,549 & Patent(s) Pending

HyFin[®] Vent Chest Seal Twin Pack

Item #	NSN#
10-0037	6510-01-642-6212



*Meets or exceeds
the current EMS
Standard of Care
and TCCC &
TECC Guidelines
for treatment of
penetrating injuries
to the chest...*



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS

— WITH —

MISSION[®]

FCP-01, 4" Emergency Bandage®

aka The Israeli Bandage®

- Best-in-class multiple use bandage addressing various wounds and wraps
- Coverage for small, medium, and large wounds
- Wrap severe wounds and small area amputations
- Can be used as an improvised tourniquet
- 4" non-adherent pad
- 8-year shelf life



Full Coverage and Multiple Functions in One Simple Package

The FCP-01 Emergency Bandage® provides everything needed to address small, medium, and large wounds in one compact package.

The versatile design is an effective solution for the control of severe bleeding, wound coverage, splinting, slings, head wraps, and creating an improvised tourniquet. The device is intuitive and simple to use for trained and untrained responders.





EMS

QUIKLOT® EMS 4X4 DRESSING AND EMS ROLLED GAUZE



QuikClot® EMS Dressings

QuikClot EMS Dressings are soft, white, nonwoven, hydrophilic gauzes impregnated with kaolin, an inorganic mineral that is both safe and effective in accelerating the body's natural clotting cascade with any exothermic reactions or use of animal or human proteins.

QuikClot® EMS 4x4 Dressing

Each QuikClot EMS 4x4 Dressing, Item #473 is packaged in a an easy-tear pouch and packed 3 per clear poly bag.

QuikClot® EMS Rolled Gauze

Each QuikClot EMS Rolled Gauze, Item #475 is packaged in an easy-tear pouch and sold individually.

QuikClot Benefits



PROVEN RESULTS

In numerous independent studies, QuikClot products have demonstrated improved time to hemostasis^{1,2} and ability to maintain a robust clot during movement.^{1,3,4}



WORKS FAST

Promotes clotting within minutes.^{1,2,5-7}



COST EFFECTIVE

Less expensive than protein-based products, with a rapid effectiveness that may reduce the need for more expensive treatments.⁸



SAFE

There are no exothermic reactions, no human or animal proteins, no thrombin, fibrinogen or shellfish products.



EASY TO USE

Intuitive to use, like standard gauze.¹ Conforms readily to the wound site and will not break down or fall apart under pressure.



For additional information, visit us at QuikClot.com



MADE IN THE USA

First Aid Pouch made of **high strength** military standard 600D nylon fabric, **scratch resistant**, **wear-resistant** and not easy to fade. Our original D ring design on the back can be used to give this IFAK bag attachment flexibility. Use them to attach a shoulder strap, attach the pouch in your car or truck. Our first aid pouch has a **2-way zipper** with silent cord pulls and it will open all the way flat so it won't snap shut while you are trying to access first aid supplies. MOLLE straps with strong snaps allow you to attach this medical pouch to any MOLLE compatible gear.



MS MARINES-THERMAL BLANKET




Also known as space blankets and emergency blankets, thermal blankets are specially crafted to trap body heat, keeping you insulated naturally. The right material, construction, and weave can easily put a thermal blanket in the same league as multiple layers of linen, insulation-wise

K K JOURA SHEAR SCISSOR



Trauma shears, also known as Tuff Cuts (a brand name), are a type of scissors used by paramedics and other emergency medical personnel to quickly and safely cut clothing from injured people.

TRAUMA CARD



Using the cards supports the key components of trauma recovery, namely becoming aware of what you are feeling, naming and symbolizing mental states with words and images, telling the story of your trauma, and understanding the meaning you took from it.

MAXSTAR MEDICAL OROPHARYNGEAL TUBE/AIRWAY.



Oropharyngeal airway devices are often used as “bite blocks” after a patient's trachea has been intubated to prevent the clenching of the teeth on the endotracheal tube. This maneuver may, however, be hazardous in children between 5 and 10 years of age with loose deciduous teeth.



Looking forward to a healthier business relationship

Thank you

SHIVANSH MEDICAL SERVICES

Plot No. G-09, Vardhman Market, Sector 2, Dwarka, New Delhi - 110075, India

Mr. Siddhant Sinha (Proprietor)

Mobile : +919350604227, +918920692107

shivanshmedicalsolutions@gmail.com

<https://shivanshsolutions.in>